

Distribution Scheme

IgG₄ Screen Nutritional 88 ELISA (Cat. No. ILE-SCG80)

88 food antigens + 6 standards + 2 controls (for 1 patient)

V11/2022

| | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 |
|----------|--------------------------|---------|--------------|-----------|--------------------|----------------------|-------------|-------------|------------------|-----------------|------------------|-----------------|
| A | 0.35 U/mL standard | lamb | sardine | camembert | plum | banana | cucumber | garlic | lentil | parsley | tea black | mustard seed |
| B | 0.7 U/mL standard | beef | plaice | parmesan | grape white&red | date | bell pepper | horseradish | soy bean | peppermint | baker's yeast | sesame |
| C | 3.5 U/mL standard | pork | tuna | roquefort | strawberry | kiwi | asparagus | red radish | mushroom | caraway | buckwheat | peanut |
| D | 17.5 U/mL standard | chicken | prawn | egg yolk | raspberry | honeydew melon | spinach | celery | cauliflower | pepper black | millet | hazelnut |
| E | 50 U/mL standard | turkey | casein | egg white | grapefruit | mango | tomato | onion | cabbage white | thyme | corn (maize) | coconut |
| F | 100 U/mL standard | cod | cow's milk | apple | orange | peach | courgette | string bean | basil | cinnamon | rice | almond |
| G | high positive control | trout | sheep's milk | pear | lemon | olive green&black | carrot | kidney bean | curry | coffee | rye | pistachio |
| H | low positive control | salmon | goat's milk | cherry | pineapple | fennel | potato | pea | dill | cocoa | wheat | walnut |