

Distribution Scheme

IgG₄ Screen Nutritional 24 ELISA (Cat. No. ILE-SCG210)

24 food antigens + 6 standards + 2 controls (for 3 patients)

V01/2023

A	0.35 U/mL standard	egg white	banana	hazelnut	0.35 U/mL standard	egg white	banana	hazelnut	0.35 U/mL standard	egg white	banana	hazelnut
B	0.7 U/mL standard	cow's milk	kiwi	peanut	0.7 U/mL standard	cow's milk	kiwi	peanut	0.7 U/mL standard	cow's milk	kiwi	peanut
C	3.5 U/mL standard	cod	strawberry	curry	3.5 U/mL standard	cod	strawberry	curry	3.5 U/mL standard	cod	strawberry	curry
D	17.5 U/mL standard	wheat	celery	pepper	17.5 U/mL standard	wheat	celery	pepper	17.5 U/mL standard	wheat	celery	pepper
E	50 U/mL standard	rye	soy bean	sesame	50 U/mL standard	rye	soy bean	sesame	50 U/mL standard	rye	soy bean	sesame
F	100 U/mL standard	barley	carrot	baker's yeast	100 U/mL standard	barley	carrot	baker's yeast	100 U/mL standard	barley	carrot	baker's yeast
G	high positive control	rice	tomato	pork	high positive control	rice	tomato	pork	high positive control	rice	tomato	pork
H	low positive control	orange	garlic	beef	low positive control	orange	garlic	beef	low positive control	orange	garlic	beef