

Distribution Scheme

IgG₄ Screen Nutritional 20 ELISA; Cat. No. ILE-SCG06

20 food antigens + 4 standards (for 4 patients)

V01/2023

A	0.35 U/mL standard	rye flour	carrot	0.35 U/mL standard	rye flour	carrot	0.35 U/mL standard	rye flour	carrot	0.35 U/mL standard	rye flour	carrot
B	0.70 U/mL standard	barley flour	tomato	0.70 U/mL standard	barley flour	tomato	0.70 U/mL standard	barley flour	tomato	0.70 U/mL standard	barley flour	tomato
C	3.5 U/mL standard	orange	hazelnut	3.5 U/mL standard	orange	hazelnut	3.5 U/mL standard	orange	hazelnut	3.5 U/mL standard	orange	hazelnut
D	17.5 U/mL standard	banana	peanut	17.5 U/mL standard	banana	peanut	17.5 U/mL standard	banana	peanut	17.5 U/mL standard	banana	peanut
E	egg white	kiwi	curry	egg white	kiwi	curry	egg white	kiwi	curry	egg white	kiwi	curry
F	cow's milk	strawberry	pepper	cow's milk	strawberry	pepper	cow's milk	strawberry	pepper	cow's milk	strawberry	pepper
G	cod	celery	sesame	cod	celery	sesame	cod	celery	sesame	cod	celery	sesame
H	wheat flour	soy bean	pork (meat)	wheat flour	soy bean	pork (meat)	wheat flour	soy bean	pork (meat)	wheat flour	soy bean	pork (meat)