

Distribution Scheme

IgG Screen Nutritional 88 ELISA (Cat. No. ILE-SCG26)

88 food antigens + 6 standards + 2 controls (for 1 patient)

V01/2023

		1	2	3	4	5	6	7	8	9	10	11
A	0.35 U/mL standard	lamb	sardine	camembert	plum	banana	cucumber	garlic	lentil	parsley	tea black	mustard seed
B	0.7 U/mL standard	beef	plaice	parmesan	grape white&red	date	bell pepper	horseradish	soy bean	peppermint	baker's yeast	sesame
C	3.5 U/mL standard	pork	tuna	roquefort	strawberry	kiwi	asparagus	red radish	mushroom	caraway	buckwheat	peanut
D	17.5 U/mL standard	chicken	prawn	egg yolk	raspberry	honeydew melon	spinach	celery	cauliflower	pepper black	millet	hazelnut
E	50 U/mL standard	turkey	casein	egg white	grapefruit	mango	tomato	onion	cabbage white	thyme	corn (maize)	coconut
F	100 U/mL standard	cod	cow's milk	apple	orange	peach	courgette	string bean	basil	cinnamon	rice	almond
G	high positive control	trout	sheep's milk	pear	lemon	olive green&black	carrot	kidney bean	curry	coffee	rye	pistachio
H	low positive control	salmon	goat's milk	cherry	pineapple	fennel	potato	pea	dill	cocoa	wheat	walnut